Dr Saurabh

JAYPEE UNIVERSITY OF INFORMATION TECHNOLOGY, WAKNAGHAT TEST -3 EXAMINATION- 2021

B.Tech V Semester (Backlog)

COURSE CODE: 15B11BT511

MAX. MARKS: 35

COURSE NAME: Bioprocess Engineering

COURSE CREDITS: 04

MAX. TIME: 2 Hours

Note: All questions are compulsory. Carrying of mobile phone during examinations will be treated as case of unfair means. For calculations, you may use simple/scientific calculator.

1. a) Determine the doubling time of a bacterial culture if it has specific growth rate 0.9 hr⁻¹. [CO3]

[1]

b) If a culture has a residence time of 12 h in a chemostat, what will be the dilution rate? [CO3]

[2]

- 2. Determine the dead band value of 5% if the pH set value of a fermenter is 8.0. [CO5] [2]
- 3. For an aerobic culture which bioreactor will you use: Fermenter A which has K_La value = 0.2 sec⁻¹ or Fermenter B having K_La value = 0.4 sec⁻¹? Give a proper justification of your selection. [CO2, CO4] [2]
- 4. An industrial fermentor containing 10,000 L of medium needs to be sterilized. The initial spore concentration in the medium is 10⁶ spores mL⁻¹. The desired probability of contamination after sterilization is 10⁻³. The death rate of spores at 121 °C is 5 min⁻¹. Assume that there is no cell death during heating and cooling phases. Determine the holding time of the sterilization process in min. [CO3, CO4] [2]
- 5. What are advantages and disadvantages of immobilization of cells? [CO1] [3]
- 6. List the major functions of the following in a fermenter: [CO5] [4]
 - i) Agitator

- ii) Baffles
- 7. Differentiate between (At least 3 main differences)

[CO1, CO5] [9]

- a) Radial flow and Axial flow Impellers
- b) Fixed and fluidized Bed Bioreactor
- c) Growth associate and Non-growth associated products

- 8. A fermentation broth with viscosity 10⁻² Pa s and density 1000 kg m⁻³ is agitated in a 2.7 m³ baffled tank using a Rushton turbine with diameter 0.5 m and stirrer speed 1 s⁻¹. Estimate the mixing time. [CO4] [3]
- 9. a) If the height-to-diameter ratio remains constant, then what happened to the surface-to-volume ratio during scaling down? [CO6] [2]
 - b) Consider the scale-up of a fermentation from a 10 l to 10,000 l vessel. The small fermenter has a height-to-diameter ratio of 3. The impeller diameter is 30% of the tank diameter. Agitator speed is 500 rpm and three Rushton impellers are used. Determine the dimensions of the large fermenter and agitator speed for Constant P/V. [CO6] [3+2]